



Special Hands
Need
Special Gloves

Baby is in safe hands.
Are you safe ?



The strength of safety
is in length of gloves

Elbowlength Gynaecology
Procedure Gloves

NULIFE and **YOU**
Together let's save lives

Introducing

Powder-free Elbow Length Gloves



Blood contact with the skin and mucous membranes of doctors occurs in 25% of vaginal deliveries and 35% of cesarean sections. **Elbow length gloves** protect the doctors forearm from significant blood and amniotic fluid contamination.

We suggest use of Powder Free gloves as:

The glove powder can enter the patients' body during surgery, which may trigger an inflammatory response by the immune system, leading to the formation of fibrous hands and post-operative adhesions.

Uterine and fallopian tube adhesions, resulting from glove powder, are a significant risk to female fertility.

Analysis of adhesions occurrence showed that 69-93% were due to foreign micro-bodies, which includes starch powder.

Powdered gloves pose to be a risk factor for post-operative wounds infections.

Available in
Powdered / Powder Free
in Sterile & non Sterile
pack



- Nitrile Examination Gloves (Non Latex)



- Examination Gloves (Latex)



- Surgical Gloves (Latex)

| Sizes | Length (mm) | Colour | Location of Thickness | | Special Features | Packing |
|--------------------------|---|---------------------|-----------------------------------|-------------|---|--|
| | | | Length Measurement | Double wall | | |
| Small Medium Large | min 350 (14) 400 (16) 450 (18) 500 (20) | Off-white to yellow | Finger (at 15mm from extreme tip) | Min. 0.45 | Powdered/ Powder free, Sterile & Non Sterile | 1 pair (1 left and 1 right) of gloves per inner wrap 1 inner wrap per pouch, 25 pouches per dispenser, 6 dispensers per carton |
| | | | Palm (at center of palm) | Min. 0.40 | | |



Manufactured & Marketed by:

MRK HEALTHCARE

B4/B5, Byculla Service Industries,
D. K. Marg, Byculla, Mumbai 400 027. INDIA
Tel: 91-22-23748371/2/3. Fax: 91-22-23748370.
E-mail: info@mrkhealthcare.com
Website: www.mrkhealthcare.com

