



Proper design and placement of drainage eyes assure maximum flow and prevents bucking of tip during insertion of catheter

Inflation lumen is built-in to avoid interior or exterior protrusion

FOLEY BALOON CATHETER

Silicon coated
natural latex catheter
for total comfort

Catheters are coated with special silicon protection to assist in smooth introduction and prevents encrustations

The taperless shaft is extremely flexible and bends easily in any direction to adjust to the patient's anatomy

It has a unique new plastic non-return valve that makes syringe insertion easy and trouble-free



NULIFE and **YOU**
Together let's save lives



Nulife Foley Catheters have been carefully designed to ensure that the balloon materials... tips and fluid channel provide the best possible performance and the highest degree of comfort to patients while in use.

Quality Assurance:

- Each foley catheter individually undergoes strict control tests before final sterilization and packaging. All balloons are inflated to ensure correct symmetry and material integrity. Additionally, every Nulife catheter is checked for surface smoothness to guard against material imperfection.
- Nulife uses high quality medical grade materials to ensure bio-compatibility.
- Special precautions are taken to ensure balloon blockage and easy, trauma-free de-catheterization.
- The balloon is made of soft rubber material for better patient comfort. The symmetrical shape, when inflated, helps to maintain excellent retention of the catheter in the bladder.

Unmatched Performance:

- Catheter balloons have been individually inflated and checked for material integrity, symmetry and integrity of the inflation/deflation system to ensure perfect functioning.
- Silicon coating gives an unmatched smooth finish for trauma-free smooth catheterization and high degree of patient comfort.
- The new plastic non-return valve makes syringe insertion easy and trouble-free even with one hand.
- Extremely flexible, easy-to-bend in any direction, the taperless shaft conforms to the patient's anatomy. Constructed with large circular drainage lumen, it achieves maximum and effective urine flow.
- Inflation lumen is built-in to avoid interior or exterior protrusion thereby, allowing the catheter to bend in any direction.
- Elliptical side eyes feature a special beveled edge. Not only does this improve smoothness and reduces trauma, but it also allows maximum fluid flow.
- Accurate placement of drainage eyes (shaped in proportion to catheter) prevent bucking of tip, making it easy to insert.

Modern Packaging :

- Nulife catheter has distinct marking on its drainage funnel, while its French size and the balloon capacity facilitates easy identification.
- Nulife Catheters have colour coded valve sleeves.
- The inner bag has been specially coated for easy opening off for

| Balloon Capacity | Sizes Available | | Packing |
|------------------|----------------------------|--------------------|-----------------|
| | Two Way | Three Way | |
| 3/5 cc | 8, 10 | - | 10 Pcs. Per Box |
| 5/15 cc | 12, 14, 16, 18, 20, 22, 24 | - | |
| 30/50 cc | 12, 14, 16, 18, 20, 22, 24 | 16, 18, 20, 22, 24 | |

Packaging: Double packed in polyethylene inner bag and outer clear view peel-open / tear-open pouch of paper and such 10 pcs. are packed in



MRK HEALTHCARE

B4/B5, Byculla Service Industries, D. K. Marg,
Mumbai 400027. Tel: 91-22-23748371/2/3.
Fax: 91-22-23748370.
E-mail: mrkhealthcare@vsnl.com
Website: www.mrkhealthcare.com

