

COMPRESSOR • NEBULIZER FOR
Aerosoltherapy
NUNEB



NULIFE and **YOU**
Together let's save lives

Breathe easy with the NUNEB range nebulizers from NULIFE

Breathe easy. With NUNEB nebulizers ,an effective treatment of respiratory diseases can be achieved.

We have a wide range of suitable devices & accessories for every individual need . Every NUNEB nebulizer is made for long lasting use while maintaining the highest standard of quality. A high aerosol output and well-controlled nebulization of precious medication make our products popular amongst many private and professional users.

NUNEB Basic

Powerful & reliable membrane type compressor specially designed for home use..
Light and compact.
Equipped with safety fuse for safety.
Safety fuse can be replaced without disassembling the instrument.

Type of compressor	Membrane
Pressure	Max 1.3 bar
Air Flow	8-10 LPM
Noise at 1 meter	62 dBA
Operating times	20 min ON, 40 min OFF
Particle sizes	0.5-5 microns
Size	344mm x 189mm x 120mm
Power Supply	(215-240vAc; 50-60 Hz)(110-130v Ac, 50-60 Hz)
Power Consumption	80 watts .

NUNEB Pro

Very powerful & reliable piston compressor nebulizer especially for professional & intensive use. Powerful compressor and matched nebulizer kit, ensures medication is atomized into fine particles faster , reaching respiratory tracts more effectively and resulting in a better medication compliance.

Equipped with safety fuse for safety.
Safety fuse can be replaced without disassembling the instrument.
Anti-shock case with handle for easy carrying.
Ample storage space for accessories.

Type of compressor	Piston
Pressure	Max 2.5 bar
Air Flow	Min 12 LPM
Noise at 1 meter	50 dBA
Operating times	60 min ON 60 min OFF
Particle sizes	0.5-5 microns
Size	344mm x 189mm x 120mm
Power Supply	(215-240vAc; 50-60 Hz)(110-130v Ac, 50-60 Hz)
Power Consumption	160 watt maximum
Weight	approx 1.9 kg

Accessories

Nebulizer kit: A perfect match for **NUNEB Basic** & **NUNEB Pro** compressors, produces particle sizes suitable for all respiratory tracts. Made from high grade boilable materials for maximum life.
Large 8ml. medication cup

- Nosepiece - Available in child & adult designs
- Mouthpiece - Available in child & adult designs
- Mask - Available in child & adult designs
- Tee - Unique design for attaching nebulizer kit with adapter with built in fresh air controller.

Substr. to change clear non toxic tubing for connecting the nebulizer cup to the compressor

MRK HEALTHCARE

B4/B5, Byculla Service Industries,
D. K. Marg, Byculla, Mumbai 400 027. INDIA
Tel: 91-22-23748371/2/3. Fax: 91-22-23748370.
E-mail: mrkhealthcare@vsnl.com
Website: www.mrkhealthcare.com

