



# **D** X-RAY OPAQUE PENROSE **R**AINAGE TUBING



**NULIFE** and **YOU**  
Together let's save lives

# X-RAY OPAQUE PENROSE DRAINAGE TUBING



## Features :

- Natural Latex Rubber material
- Smooth tubing wall both inside & on the outside
- Sterile ready to implant
- Complies with ASTM D 3579
- Sizes available to suit various requirements

| Item                                 | Available In Length | Sizes       | Box of Qty. |
|--------------------------------------|---------------------|-------------|-------------|
| X-Ray Opaque Penrose Drainage Tubing | 18"/12"             | 1/4" , 3/8" | 50 Pcs.     |
|                                      |                     | 1/2" , 5/8" |             |
|                                      |                     | 3/4" , 7/8" |             |
|                                      |                     | 1" , 1 1/4" |             |

### Packaging

- Every piece of penrose drain is folded & packed in a peel off pouch.
- 50 such pouches are packed in a dispenser box.
- Sterilization by Ethylene Oxide.



## MRK HEALTHCARE

B4/B5, Byculla Service Industries, D. K. Marg,  
Mumbai 400 027. Tel: 91-22-23748371/2/3.  
Fax: 91-22-23748370.  
E-mail: mrkhealthcare@vsnl.com  
Website: www.mrkhealthcare.com

