



NULIFE
Since 1965

NULIFE THREWAY STOPCOCK & EXTENSION TUBING



EXTENSION TUBING



EXTENSION TUBING WITH
INTEGRATED THREE WAY STOP COCK



HIGH PRESSURE
EXTENSION TUBING



Patient care devices
with excellence



NULIFE and **YOU**
Together let's save lives



THREWAY STOPCOCK - NULIFE®

- Accurate drug administration.
- Can withstand the pressure in normal procedures
- Arrow indication marks on the tap indicates direction of flow.
- Leak proof up to three bars.
- Rotating luer lock and threaded female ports facilitate safe & secure connection
- Internationally accepted 6% luer taper
- Lipid resistant is also available.
- Disposable, Sterile and Non-

EXTENSION TUBING - NULIFE®

- Extends the pathway between the infusion site and the source of infusion.
- Helps to take manipulation site away from the insertion site.
- Male - Male connection at both ends is also available.
- Considerably reduces the chances of mechanical irritation and infection thereby reducing trauma to the patient.
- Made from non-toxic medical grade PVC material.
- Disposable, Sterile and Non-Pyrogenic.

EXTENSION TUBING with THREWAY STOPCOCK NULIFE®

- Assist in the administration of intravenous infusion.
- Extends the pathway between the infusion site and the source of infusion.
- Facilitates simultaneous administration of additional fluid from side port.
- Helps to take manipulation site away from the insertion site.
- Various lengths offer flexibility of use according to applications.
- Fitted with male luer lock at one end & the three way stop cock at the other.
- Disposable, Sterile and Non-

HIGH PRESSURE EXTENSION TUBING - NULIFE®

- Fitted with male & female luer locks for sure connection.
- Does not bulge even at high pressure.
- Minimal dead space.
- Withstand a pressure of 40 kg/cm².
- Disposable, Sterile and Non-Pyrogenic.

PACKING DETAILS

THREWAY STOPCOCK			
CODE (Normal)	CODE (Lipid resistant)	Tap Colour	Packing
10012 10022 10032 10042	20012 20022 20032 20042	Blue White Red White	50 / 500

EXTENSION TUBING							
Pouch pack				Soft blister pack			
CODE	Ext.-Int. Dia.	Length Cm	Packing	CODE	Ext.-Int. Dia.	Length Cm	Packing
60051 60071 60101 60251 60501 60751 61001 61251	4.0-2.8 4.0-2.8 4.0-2.8 4.0-2.8 4.0-2.8 4.0-2.8 4.0-2.8 4.0-2.8	5 7 10 25 50 75 100 125	50/1000	60053 60073 60103 60253 60503 60753 61003 61253	4.0-2.8 4.0-2.8 4.0-2.8 4.0-2.8 4.0-2.8 4.0-2.8 4.0-2.8 4.0-2.8	5 7 10 25 50 75 100 125	25/500
61501 62001	4.0-2.8 4.0-2.8	150 200	50/750	61503 62003	4.0-2.8 4.0-2.8	150 200	25/300

EXTENSION TUBING with THREWAY STOPCOCK							
Pouch pack				Soft blister pack			
CODE	Ext.-Int. Dia.	Length Cm	Packing	CODE	Ext.-Int. Dia.	Length Cm	Packing
70051 70071 70101 70251 70501 70751	4.5-2.8 4.5-2.8 4.5-2.8 4.5-2.8 4.5-2.8 4.5-2.8	5 7 10 25 50 75	25/1000 25/750	70053 70073 70103 70253 70503 70753	4.5-2.8 4.5-2.8 4.5-2.8 4.5-2.8 4.5-2.8 4.5-2.8	5 7 10 25 50 75	25/500
71001 71251 71501 72001	4.5-2.8 4.5-2.8 4.5-2.8 4.5-2.8	100 125 150 200	25/500	71003 71253 71503 72003	4.5-2.8 4.5-2.8 4.5-2.8 4.5-2.8	100 125 150 200	25/300

HIGH PRESSURE EXTENSION TUBING							
Pouch pack				Soft blister pack			
CODE	Ext.-Int. Dia.	Length Cm	Packing	CODE	Ext.-Int. Dia.	Length Cm	Packing
90051 90071 90101 90251 90501 90751 91001 91251	3.0/1.0 3.0/1.0 3.0/1.0 3.0/1.0 3.0/1.0 3.0/1.0 3.0/1.0 3.0/1.0	5 7 10 25 50 75 100 125	50/1000	90053 90073 90103 90253 90503 90753 91003 91253	3.0/1.0 3.0/1.0 3.0/1.0 3.0/1.0 3.0/1.0 3.0/1.0 3.0/1.0 3.0/1.0	5 7 10 25 50 75 100 125	25/500
91501 92001	3.0/1.0 3.0/1.0	150 200	50/750	91503 92003	3.0/1.0 3.0/1.0	150 200	25/300



MRK HEALTHCARE

B4/B5, Byculla Service Industries, D. K. Marg,
Mumbai 400 027. Tel: 91-22-23748371/2/3.
Fax: 91-22-23748370.
E-mail: mrkhealthcare@vsnl.com

