



Snug fit without discomfort or constriction

Individual elastic foam tape is provided to hold catheter and to prevent slipping

One piece catheter

U DRAIN

The comfortable, easy-to-use Male Incontinence Device...

...with specially compounded rubber for user convenience.

Reinforced funnel prevents twisting and collapse

Kink resistant bulb maintains free flow of urine and prevents its retention in catheter

Long connector with large diameter provides maximum drainage and device fits into any standard urine bag or receptor



NULIFE and **YOU**
Together let's save lives

U-DRAIN



Nulife U-Drain is a male incontinence device that serves as a partial alternative to catheterization. It can be used as a replacement for catheters in many cases, by males who have lost control over urination or by those who need to urinate frequently. It can also be used by bed-ridden patients.

Its simple and easy-to-use design ensures that patients can use

Nulife U-Drain is very comfortable to use, being as soft as a hat, made of soft, comfortable rubber that provides a snug fit without constriction. It is flexible enough to accommodate changes in penis size during use. Nulife provides an individual adhesive strip with



Item	Size	Diameter	Packing
U-DRAIN (Male incontinence Device)	Extra Small	20 mm	50 Pcs Per Box
	Small	25 mm	50 Pcs Per Box
	Medium	30 mm	50 Pcs Per Box
	Large	35 mm	50 Pcs Per Box
	Extra Large	40 mm	50 Pcs Per Box

Packing: Each piece is packed along with adhesive strip in a peel-off / tear-off paper pouch. A box contains 50 such pcs.



MRK HEALTHCARE

B4/B5, Byculla Service Industries, D. K. Marg,
Mumbai 400027. Tel: 91-22-23748371/2/3.
Fax: 91-22-23748370.
E-mail: mrkhealthcare@vsnl.com
Website: www.mrkhealthcare.com

